



RNF26 Antibody, HRP conjugated

Product Code	CSB-PA019883LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY78
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF26 protein (246-433AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF26 (EC 2.3.2.27) (RING finger protein 26), RNF26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RNF26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.