



RNF31 Antibody, Biotin conjugated

Product Code	CSB-PA856927LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EP0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF31 protein (421-700AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF31 (EC 2.3.2.27) (HOIL-1-interacting protein) (HOIP) (RING finger protein 31) (RING-type E3 ubiquitin transferase RNF31) (Zinc in-between-RING-finger ubiquitin-associated domain protein), RNF31, ZIBRA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.