




MRM3 Antibody

Product Code	CSB-PA867196LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC36
Immunogen	Recombinant Human rRNA methyltransferase 3, mitochondrial protein (41-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	rRNA methyltransferase 3, mitochondrial (EC 2.1.1.-) (16S rRNA (guanosine(1370)-2'-O)-methyltransferase) (16S rRNA [Gm1370] 2'-O-methyltransferase) (RNA methyltransferase-like protein 1), MRM3, RNMTL1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRM3

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA867196LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.