



ROGDI Antibody, Biotin conjugated

Product Code	CSB-PA867108LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZN7
Immunogen	Recombinant Human Protein rogdi homolog protein (14-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein rogdi homolog, ROGDI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ROGDI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.