



RPL10L Antibody, FITC conjugated

Product Code	CSB-PA850330LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96L21
Immunogen	Recombinant Human 60S ribosomal protein L10-like protein (1-214AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L10-like (Large ribosomal subunit protein uL16-like), RPL10L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL10L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.