



RPRD1B Antibody, Biotin conjugated

Product Code	CSB-PA878854HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQG5
Immunogen	Recombinant Human Regulation of nuclear pre-mRNA domain-containing protein 1B protein (2-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulation of nuclear pre-mRNA domain-containing protein 1B (Cell cycle-related and expression-elevated protein in tumor), RPRD1B, C20orf77 CREPT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RPRD1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.