



RPRD1B Antibody, HRP conjugated

Product Code	CSB-PA878854HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQG5
Immunogen	Recombinant Human Regulation of nuclear pre-mRNA domain-containing protein 1B protein (2-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulation of nuclear pre-mRNA domain-containing protein 1B (Cell cycle-related and expression-elevated protein in tumor), RPRD1B, C20orf77 CREPT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RPRD1B