





RPUSD4 Antibody

Product Code	CSB-PA836192LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CM3
Immunogen	Recombinant Human Mitochondrial RNA pseudouridine synthase RPUSD4 protein (41-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial RNA pseudouridine synthase RPUSD4 (EC 5.4.99.-) (RNA pseudouridylate synthase domain-containing protein 4), RPUSD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPUSD4


Image



Western blot
 All lanes: RPUSD4 antibody at 0.5µg/ml +
 MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 43, 39 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA836192LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA836192LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.