



RPUSD4 Antibody, Biotin conjugated

Product Code	CSB-PA836192LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CM3
Immunogen	Recombinant Human Mitochondrial RNA pseudouridine synthase RPUSD4 protein (41-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial RNA pseudouridine synthase RPUSD4 (EC 5.4.99.-) (RNA pseudouridylate synthase domain-containing protein 4), RPUSD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPUSD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.