



CNOT9 Antibody, FITC conjugated

Product Code	CSB-PA842140LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92600
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 9 protein (1-258AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCR4-NOT transcription complex subunit 9 (Cell differentiation protein RQCD1 homolog) (Rcd-1), CNOT9, RCD1 RQCD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CNOT9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.