





RSPH14 Antibody

Product Code	CSB-PA887042LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHP6
Immunogen	Recombinant Human Radial spoke head 14 homolog protein (1-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radial spoke head 14 homolog (Rhabdoid tumor deletion region protein 1), RSPH14, RTDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RSPH14

Image



Western blot
 All lanes: RSPH14 antibody at 8µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 39 kDa
 Observed band size: 39 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA887042LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.