



RSPH14 Antibody, HRP conjugated

Product Code	CSB-PA887042LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHP6
Immunogen	Recombinant Human Radial spoke head 14 homolog protein (1-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radial spoke head 14 homolog (Rhabdoid tumor deletion region protein 1), RSPH14, RTDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RSPH14