



RTBDN Antibody

Product Code	CSB-PA020560LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSG5
Immunogen	Recombinant Human Retbindin protein (63-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retbindin, RTBDN
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTBDN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.