



RTCA Antibody, Biotin conjugated

Product Code	CSB-PA520755LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00442
Immunogen	Recombinant Human RNA 3'-terminal phosphate cyclase protein (167-379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA 3'-terminal phosphate cyclase (RNA cyclase) (RNA-3'-phosphate cyclase) (EC 6.5.1.4) (RNA terminal phosphate cyclase domain-containing protein 1) (RTC domain-containing protein 1), RTCA, RPC RPC1 RTC1 RTCD1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTCA