



RTCA Antibody, FITC conjugated

Product Code	CSB-PA520755LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00442
Immunogen	Recombinant Human RNA 3'-terminal phosphate cyclase protein (167-379AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA 3'-terminal phosphate cyclase (RNA cyclase) (RNA-3'-phosphate cyclase) (EC 6.5.1.4) (RNA terminal phosphate cyclase domain-containing protein 1) (RTC domain-containing protein 1), RTCA, RPC RPC1 RTC1 RTCD1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.