





RTN4IP1 Antibody

Product Code	CSB-PA848832LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWV3
Immunogen	Recombinant Human Reticulon-4-interacting protein 1, mitochondrial protein (41-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4-interacting protein 1, mitochondrial (NOGO-interacting mitochondrial protein), RTN4IP1, NIMP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTN4IP1

Image



Immunohistochemistry of paraffin-embedded human appendix tissue using CSB-PA848832LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA848832LA01HU at dilution of 1:100