



RUNDC3A Antibody, FITC conjugated

Product Code	CSB-PA685437LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59EK9
Immunogen	Recombinant Human RUN domain-containing protein 3A protein (1-405AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RUN domain-containing protein 3A (Rap2-interacting protein 8) (RPIP-8), RUNDC3A, RAP2IP RPIP8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RUNDC3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.