



RUSC1 Antibody, FITC conjugated

Product Code	CSB-PA871585LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVN2
Immunogen	Recombinant Human RUN and SH3 domain-containing protein 1 protein (231-433AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RUN and SH3 domain-containing protein 1 (New molecule containing SH3 at the carboxy-terminus) (Nesca), RUSC1, NESCA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RUSC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.