



S100PBP Antibody, Biotin conjugated

Product Code	CSB-PA822185LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BU1
Immunogen	Recombinant Human S100P-binding protein (276-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S100P-binding protein (S100P-binding protein Riken), S100PBP, S100PBPR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S100PBP