



SAMD3 Antibody, HRP conjugated

Product Code	CSB-PA843281LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6K7
Immunogen	Recombinant Human Sterile alpha motif domain-containing protein 3 protein (62-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sterile alpha motif domain-containing protein 3 (SAM domain-containing protein 3), SAMD3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SAMD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.