



SAMSN1 Antibody, Biotin conjugated

Product Code	CSB-PA885722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSI8
Immunogen	Recombinant Human SAM domain-containing protein SAMSN-1 protein (1-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SAM domain-containing protein SAMSN-1 (Hematopoietic adaptor containing SH3 and SAM domains 1) (Nash1) (SAM domain, SH3 domain and nuclear localization signals protein 1) (SH3-SAM adaptor protein), SAMSN1, HACS1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SAMSN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.