



SART3 Antibody, HRP conjugated

Product Code	CSB-PA624010HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15020
Immunogen	Recombinant Human Squamous cell carcinoma antigen recognized by T-cells 3 protein (600-900AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Squamous cell carcinoma antigen recognized by T-cells 3 (SART-3) (Tat-interacting protein of 110 kDa) (Tip110) (p110 nuclear RNA-binding protein), SART3, KIAA0156 TIP110
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SART3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.