



SCAMP5 Antibody, HRP conjugated

Product Code	CSB-PA855032LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAC9
Immunogen	Recombinant Human Secretory carrier-associated membrane protein 5 protein (170-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretory carrier-associated membrane protein 5 (Secretory carrier membrane protein 5) (hSCAMP5), SCAMP5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCAMP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.