



SCAND2P Antibody, Biotin conjugated

Product Code	CSB-PA861977LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZW5
Immunogen	Recombinant Human Putative SCAN domain-containing protein SCAND2P protein (1-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative SCAN domain-containing protein SCAND2P (SCAN domain-containing protein 2 pseudogene), SCAND2P, SCAND2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCAND2P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.