



SCNN1G Antibody, Biotin conjugated

Product Code	CSB-PA020851LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51170
Immunogen	Recombinant Human Amiloride-sensitive sodium channel subunit gamma protein (90-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amiloride-sensitive sodium channel subunit gamma (Epithelial Na(+)) channel subunit gamma) (ENaCG) (Gamma-ENaC) (Gamma-NaCH) (Nonvoltage-gated sodium channel 1 subunit gamma) (SCNEG), SCNN1G
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCNN1G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.