



SCRG1 Antibody, HRP conjugated

Product Code	CSB-PA530448LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75711
Immunogen	Recombinant Human Scrapie-responsive protein 1 protein (21-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scrapie-responsive protein 1 (Scrapie-responsive gene 1 protein) (ScRG-1), SCRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCRG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.