



Scrg1 Antibody, HRP conjugated

Product Code	CSB-PA526023HB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O88745
Immunogen	Recombinant Mouse Scrapie-responsive protein 1 protein (21-98AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scrapie-responsive protein 1 (Scrapie-responsive gene 1 protein) (ScRG-1), Scrg1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Scrg1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.