



SDCBP2 Antibody, FITC conjugated

Product Code	CSB-PA872430LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H190
Immunogen	Recombinant Human Syntenin-2 protein (11-207AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntenin-2 (Similar to TACIP18) (SITAC) (Syndecan-binding protein 2), SDCBP2, SITAC18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SDCBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.