



SDHC Antibody, FITC conjugated

Product Code	CSB-PA859106LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99643
Immunogen	Recombinant Human Succinate dehydrogenase cytochrome b560 subunit, mitochondrial protein (30-62AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Succinate dehydrogenase cytochrome b560 subunit, mitochondrial (Integral membrane protein CII-3) (QPs-1) (QPs1) (Succinate dehydrogenase complex subunit C) (Succinate-ubiquinone oxidoreductase cytochrome B large subunit) (CYBL), SDHC, CYB560 SDH3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SDHC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.