



SEC22A Antibody, HRP conjugated

Product Code	CSB-PA856954LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IW7
Immunogen	Recombinant Human Vesicle-trafficking protein SEC22a protein (10-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-trafficking protein SEC22a (SEC22 vesicle-trafficking protein homolog A) (SEC22 vesicle-trafficking protein-like 2), SEC22A, SEC22L2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC22A