



SEC24C Antibody, HRP conjugated

Product Code	CSB-PA020952LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53992
Immunogen	Recombinant Human Protein transport protein Sec24C protein (200-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein transport protein Sec24C (SEC24-related protein C), SEC24C, KIAA0079
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC24C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.