



SEC31A Antibody, FITC conjugated

Product Code	CSB-PA020954LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94979
Immunogen	Recombinant Human Protein transport protein Sec31A protein (335-503AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein transport protein Sec31A (ABP125) (ABP130) (SEC31-like protein 1) (SEC31-related protein A) (Web1-like protein), SEC31A, KIAA0905 SEC31L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC31A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.