



Serpinh1 Antibody, FITC conjugated

Product Code	CSB-PA021087LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19324
Immunogen	Recombinant Mouse Serpin H1 protein (18-417AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin H1 (47 kDa heat shock protein) (Collagen-binding protein) (Colligin) (Serine protease inhibitor J6), Serpinh1, Cbp1 Hsp47
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Serpinh1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.