



Serpinh1 Antibody, FITC conjugated

Product Code	CSB-PA021087LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19324
Immunogen	Recombinant Mouse Serpin H1 protein (18-417AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin H1 (47 kDa heat shock protein) (Collagen-binding protein) (Colligin) (Serine protease inhibitor J6), Serpinh1, Cbp1 Hsp47
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Serpinh1