



SERTAD1 Antibody, Biotin conjugated

Product Code	CSB-PA883406LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHV2
Immunogen	Recombinant Human SERTA domain-containing protein 1 protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SERTA domain-containing protein 1 (CDK4-binding protein p34SEI1) (SEI-1) (Transcriptional regulator interacting with the PHD-bromodomain 1) (TRIP-Br1), SERTAD1, SEI1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SERTAD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.