



SETBP1 Antibody

Product Code	CSB-PA896935LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6X0
Immunogen	Recombinant Human SET-binding protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SET-binding protein (SEB), SETBP1, KIAA0437
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	SETBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.