



SETX Antibody, Biotin conjugated

Product Code	CSB-PA800097LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z333
Immunogen	Recombinant Human Probable helicase senataxin protein (2401-2677AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable helicase senataxin (EC 3.6.4.-) (Amyotrophic lateral sclerosis 4 protein) (SEN1 homolog) (Senataxin), SETX, ALS4 KIAA0625 SCAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SETX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.