




SGMS2 Antibody

Product Code	CSB-PA854159LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHU3
Immunogen	Recombinant Human Phosphatidylcholine:ceramide cholinephosphotransferase 2 protein (1-79AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGMS2

Image



Western Blot

Positive WB detected in: A549 whole cell lysate(30µg), HepG2 whole cell lysate(30µg), Hep2 whole cell lysate(30µg), PC-3 whole cell lysate(30µg), Mouse Spleen tissue lysate(30µg), Mouse Lung tissue lysate(30µg), Rat Lung tissue lysate(30µg), Rat Spleen tissue lysate(30µg)

All lanes: SGMS2 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 43 kDa

Observed band size: 43 kDa

Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.