



SGTB Antibody, Biotin conjugated

Product Code	CSB-PA846614LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EQ0
Immunogen	Recombinant Human Small glutamine-rich tetratricopeptide repeat-containing protein beta protein (14-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small glutamine-rich tetratricopeptide repeat-containing protein beta (Beta-SGT) (Small glutamine-rich protein with tetratricopeptide repeats 2), SGTB, SGT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGTB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.