



# SGTB Antibody, HRP conjugated

Product Code	CSB-PA846614LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EQ0
Immunogen	Recombinant Human Small glutamine-rich tetratricopeptide repeat-containing protein beta protein (14-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small glutamine-rich tetratricopeptide repeat-containing protein beta (Beta-SGT) (Small glutamine-rich protein with tetratricopeptide repeats 2), SGTB, SGT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGTB