



SH3BGR Antibody, Biotin conjugated

Product Code	CSB-PA021219LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55822
Immunogen	Recombinant Human SH3 domain-binding glutamic acid-rich protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding glutamic acid-rich protein (SH3BGR protein) (21-glutamic acid-rich protein) (21-GARP), SH3BGR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SH3BGR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.