



SHISA5 Antibody, Biotin conjugated

Product Code	CSB-PA818677LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N114
Immunogen	Recombinant Human Protein shisa-5 protein (56-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein shisa-5 (Putative NF-kappa-B-activating protein 120) (Scotin), SHISA5, SCOTIN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SHISA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.