



SIGLEC15 Antibody, Biotin conjugated

Product Code	CSB-PA761623LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZMC9
Immunogen	Recombinant Human Sialic acid-binding Ig-like lectin 15 protein (20-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid-binding Ig-like lectin 15 (Siglec-15) (CD33 antigen-like 3), SIGLEC15, CD33L3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SIGLEC15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.