



SIGLEC15 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA761623LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6ZMC9 |
| Immunogen | Recombinant Human Sialic acid-binding Ig-like lectin 15 protein (20-263AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sialic acid-binding Ig-like lectin 15 (Siglec-15) (CD33 antigen-like 3), SIGLEC15, CD33L3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | SIGLEC15 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |