



SIGLEC6 Antibody, Biotin conjugated

Product Code	CSB-PA021300LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43699
Immunogen	Recombinant Human Sialic acid-binding Ig-like lectin 6 protein (76-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid-binding Ig-like lectin 6 (Siglec-6) (CD33 antigen-like 1) (CDw327) (Obesity-binding protein 1) (OB-BP1) (CD antigen CD327), SIGLEC6, CD33L CD33L1 OBBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SIGLEC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.