



SIMC1 Antibody, Biotin conjugated

Product Code	CSB-PA822768LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDZ2
Immunogen	Recombinant Human SUMO-interacting motif-containing protein 1 protein (18-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SUMO-interacting motif-containing protein 1, SIMC1, C5orf25
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SIMC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.