



SLA Antibody, FITC conjugated

Product Code	CSB-PA613387LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13239
Immunogen	Recombinant Human Src-like-adapter protein (2-276AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Src-like-adapter (Src-like-adapter protein 1) (SLAP-1) (hSLAP), SLA, SLAP SLAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.