




SLC12A9 Antibody

Product Code	CSB-PA866307LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXP2
Immunogen	Recombinant Human Solute carrier family 12 member 9 protein (126-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 9 (Cation-chloride cotransporter 6) (hCCC6) (Cation-chloride cotransporter-interacting protein 1) (CCC-interacting protein 1) (hCIP1) (Potassium-chloride transporter 9) (WO3.3), SLC12A9, CCC6 CIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC12A9

Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA866307LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)