




# SLC12A9 Antibody

<b>Product Code</b>	CSB-PA866307LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BXP2
<b>Immunogen</b>	Recombinant Human Solute carrier family 12 member 9 protein (126-173AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IF; Recommended dilution: IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Solute carrier family 12 member 9 (Cation-chloride cotransporter 6) (hCCC6) (Cation-chloride cotransporter-interacting protein 1) (CCC-interacting protein 1) (hCIP1) (Potassium-chloride transporter 9) (WO3.3), SLC12A9, CCC6 CIP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	SLC12A9

## Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA866307LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.