



SLC22A18AS Antibody, Biotin conjugated

Product Code	CSB-PA847610LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1D0
Immunogen	Recombinant Human Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein (Organic cation transporter-like protein 2 antisense protein) (Solute carrier family 22 member 1-like antisense protein) (Solute carrier family 22 member 18 antisense protein) (p27-Beckwith-Wiedemann region 1 B) (p27-BWR1B), SLC22A18AS, BWR1B BWSCR1B ORCTL2S SLC22A1LS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SLC22A18AS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.